



The CXI File Format for Coherent X-ray Imaging

Version 1.6

Filipe R. N. C. Maia
December 17, 2019

Contents

1	Introduction	1
1.1	Goals	1
2	The Design of the CXI	1
2.1	HDF5	1
2.2	NeXus	2
2.3	CXI	2
3	CXI by example	3
3.1	A minimal CXI file	3
3.2	A typical raw data CXI file	3
3.3	A CXI file with raw tomographic data	5
3.4	CXI file from a ptychography experiment	6
3.5	CXI from a modular detector	7
3.6	CXI stack of images from a modular detector	8
3.7	CXI file after initial analysis of a diffraction image	9
3.8	CXI file of a phased 3D image	10
4	NeXus compatibility	11
5	Datatypes	12
5.1	Complex Numbers	12
6	Default units for CXI entries	13
6.1	Angles	13
6.2	Dates	13
7	Memory Layout	14
8	Geometry	14
8.1	Coordinate System	14
8.2	The local coordinate system of objects	14
8.2.1	The orientation of pixel detectors	15
9	Scans	16
9.1	Dimension Scales	17
9.2	Implicit Axes	17
9.3	Modular Pixel Detectors	17

A	CXI entries reference	19
A.1	Top level (root)	19
A.2	Attenuator	19
A.3	Data	20
A.4	Detector	20
A.5	Entry	23
A.6	Geometry	23
A.7	Image	24
A.8	Instrument	27
A.9	Monochromator	27
A.10	Note	28
A.11	Process	28
A.12	Result	29
A.13	Sample	31
A.14	Source	32
B	Changes from previous versions	33
B.1	Version 1.5	33
B.2	Version 1.4	33
B.3	Version 1.3	33
B.4	Version 1.2.1	34
B.5	Version 1.2	34
B.6	Version 1.1	34
B.7	Version 1.0	34
C	Diagrams color code	35
D	Code examples	36
D.1	Creating a minimal CXI file	36
D.2	Creating a typical raw CXI file	37
D.3	NeXus compatible version of the typical raw CXI file	39
D.4	3D complex valued Image file	41
D.5	Spartan viewer for CXI files	43

1 Introduction

The CXI file format was created as common format for all the data in the Coherent X-ray Imaging Data Bank (CXIDB). Naturally its scope is all experimental data collected during Coherent X-ray Imaging experiments as well as all data generated during the analysis of the experimental data.

1.1 Goals

The CXI file format aims to create a data format with the following requirements:

1. Simple - both writing and reading should be made simple.
2. Flexible - users should be able to easily extend it.
3. Fast - it should be efficient so as not to become a bottleneck.
4. Extendable - new features should be easily added without breaking compatibility with previous versions.
5. Unambiguous - it should be possible to interpret the files without using external information.
6. Compatible - the format should be as compatible as possible with existing formats

These are hard and often contradicting requirements (e.g. simple and unambiguous). When such conflicts occur the highest ranked goal is often preferred (e.g. simple).

2 The Design of the CXI

2.1 HDF5

The HDF5 format is the basis of CXI format. CXI is not really a completely new file format but simply a set of rules designed to create HDF5 files with a common structure and to allow a uniform and consistent interpretation of such files.

HDF5 was chosen as the basis because it is a widely used high performance scientific data format which many programs can already, at least partially, read and write. It also brings with it the almost automatic fulfilment of requirements 2, 3 and 4. HDF5 version 1.8 or higher is required as previous versions don't support all features required by CXI.

2.2 NeXus

Another important influence in the design of CXI is the NeXus file format for neutron, x-ray and muon science. As CXI, NeXus is also based on HDF5 (although it can use others basis formats such as XML) and shares many of the goals of CXI, with one important exception - it is *not* a simple format. All NeXus file require attributes in HDF5 files and using many existing programs it is laborious and hard, if not outright impossible, to create such attributes. The large scope of NeXus also results in files that can be complex to read and interpret.

NeXus classes are the fundamental pieces that make up NeXus files. The documentation for the classes can be found in <http://download.nexusformat.org/doc/html/ClassDefinitions.html>. Each class has a set of members which can either be data fields with well defined members (e.g. the `data` field of the `NXdata` class), or they can be other classes in which case the name is arbitrary. Each NeXus class corresponds to an HDF5 group with one or more attributes (class name, units, long name, etc...) and each class data field corresponds to an HDF5 dataset.

2.3 CXI

The approach used in CXI is to use HDF5 as basis format and adopt a small set of rules derived from the NeXus format with some additional restraints which should make the format simpler to write and interpret. This allows us to relax the format, making it simple for new users, while still retaining enough information to reconstruct an NeXus compile file if necessary.

Below are the main rules used in CXI files:

- Groups representing NeXus classes should be named like the class with the `NX` prefix removed and with `_N` added where `N` are consecutive positive integers starting at 1 (e.g. `entry_1` represents the first `NXentry`).
- NeXus class fields should be used as much as possible.
- Each data field must have default units.
- Non default units are defined using the `units` attribute of relevant dataset, just like in NeXus.

3 CXI by example

3.1 A minimal CXI file

Figure 1 shows a diagram a minimal CXI file. Each measurement (which is usually an exposure/shot) in a CXI file is stored inside a group named `entry_N` where N is a positive integer. The first such entry should be named `entry_1` the second `entry_2` and so on. Each entry is a member of the `NXentry` class in NeXus language.




Figure 1: Diagram of a minimal CXI file with a single image.

The arrangements of the HDF5 groups is similar to that in a NeXus file. The main difference is that in a NeXus file entries are not identified by name but by the `NX_class` attribute which is set to `NXentry`. In a CXI file we restrict the possible names of a measurement group to only `entry_N` simplifying things by not requiring the attribute. Yet a simple post processing program could take a CXI file and easily convert it into a NeXus file by adding the required attributes. This is the *main design idea* of the CXI format and we will see examples of it in many places.

Each entry can have one or more data groups. In this case we only had one detector so we only have one data group named `data_1`. The raw data is stored in the `data` field of the `data_1` group. As no units are specified the data is assumed to be in “counts” (see 6). Also no experimental data about the experimental conditions is stored.

3.2 A typical raw data CXI file




Figure 2: Diagram of a typical CXI file for storing raw data from a single shot.

3.3 A CXI file with raw tomographic data

This file exemplifies the use of Scans.




Figure 3: Diagram of a CXI file for storing raw tomographic data.

3.4 CXI file from a ptychography experiment




Figure 4: Diagram of a CXI file with 256 diffraction images recorded on a 1024x1024 pixels CCD detector from a ptychography experiment.

3.5 CXI from a modular detector

This file exemplifies the use of [Modular Pixel Detector](#).




Figure 5: Diagram of a CXI file with an image from a pixel detector made of 64 modules, such as the [LCLS CSPAD](#).

3.6 CXI stack of images from a modular detector

This file exemplifies the use of **Modular Pixel Detector** together with **Scans** to store a large dataset of 30,000 images.




Figure 6: Diagram of a CXI file with a stack of 30,000 images from a pixel detector made of 64 modules, such as the LCLS CSPAD.

3.7 CXI file after initial analysis of a diffraction image




Figure 7: Diagram of a CXI file of an analysed image from a single particle diffraction experiment.

3.8 CXI file of a phased 3D image




Figure 8: Diagram of a CXI file of a phased 3D image from a single particle diffraction experiment.

4 NeXus compatibility

One of the most appealing features of the CXI format is that while it is relatively simple it can still be unambiguously converted to follow the NeXus format.

Figure 9 shows the conversion of the example in figure 2 to the NeXus format. Note that the only change are the addition of several attributes that are necessary for NeXus. This means that CXI files can be easily converted to NeXus, then edited using NeXus based tools if necessary, and then read back using a CXI based tool.




Figure 9: Diagram of a NeXus compatible CXI file for storing raw data from a single shot.

5 Datatypes

HDF5 covers a large variety of native datatypes including integers, floating point numbers and character string(including UTF-8 support). It also takes care of the conversion of datatypes when reading and writing files (eliminating endianness problems for example).

Most of the data should be saved in the same format as it was created/acquired. For example CCD images acquired as 16 bit integers should be saved using the `H5T_NATIVE_SHORT` HDF5 Type.

5.1 Complex Numbers

A notable omission of the HDF5 1.8 standard, in which this format is based on, is the lack of a standard way to store complex numbers.

The CXI convention for storing complex numbers is to use a compound data type with two elements named `r` and `i`. The real part of the number should obviously be saved in the element named `r` and the imaginary part in the one named `i`. This follows the convention adopted by PyTables as well as h5py.

Below you can see a C99 code snippet showing how to create a CXI compatible HDF5 compound type for double complex numbers.

Listing 1: Creating a double complex type

```
hid_t complex_id = H5Tcreate(H5T_COMPOUND,
                            sizeof(double complex));
H5Tinsert(complex_id, "r", 0, H5T_NATIVE_DOUBLE);
H5Tinsert(complex_id, "i", sizeof(double), H5T_NATIVE_DOUBLE);
```

6 Default units for CXI entries

The default units for all CXI entries are SI base units (see table 1) with *no exceptions*.

Table 1: SI (and common derived) base units for different quantities

Quantity	Units	Abbreviation
length	meter	m
mass	kilogram	kg
time	second	s
electric current	ampere	A
temperature	kelvin	K
amount of substance	mole	mol
luminous intensity	candela	cd
frequency	hertz	Hz
force	newton	N
pressure	pascal	Pa
energy	joule	J
power	watt	W
electric potential	volt	V
capacitance	farad	F
electric resistance	ohm	Ω
absorbed dose	gray	Gy
area	square meter	m ²
volume	cubic meter	m ³

6.1 Angles

Angles are always defined in degrees *not* in radians.

6.2 Dates

Dates are always specified according to the [ISO 8601](#). This means for example “1996-07-31T21:15:22+0600”. Note the “T” separating the data from the time and the “+0600” timezone specification.

All of these are mandatory. They derive from the use of ISO 8601 in NeXus. This way compatibility is ensured.

7 Memory Layout

All multidimensional arrays must be stored with the fastest changing dimension being the last dimension, and the slowest changing dimension being the first dimension, also known as row major. This is defined in the HDF5 standard.

8 Geometry

8.1 Coordinate System

CXI uses the same coordinate system as NeXus which itself is based on the McStas coordinate system (NeXus User Manual section 2.2.1).

The CXI coordinate system is a right handed system. The z axis is parallel to the X-ray beam, with the positive z direction pointing away from the light source, in the downstream direction. This is the *opposite* of the definition of the International Tables for Crystallography, volume G. The y axis is vertical with the positive direction pointing up, while the x axis is horizontal completing the right handed system (see Fig. 10).

The origin of the CXI coordinate system is defined by the point where the X-ray beam meets the sample.




Figure 10: The coordinate system used by CXI. The intersection of the X-ray beam with the sample define the origin of the system. The z axis is parallel to the beam and points downstream.

8.2 The local coordinate system of objects

For many detectors their location and orientation is crucial to interpret results.

Translations and rotations are used to define the absolute position of each object. But to be able to apply these transformations we need to know what is the origin of the local coordinate system of each object.

Unless otherwise specified the origin should be assumed to be the geometrical center of the object in question. The default orientation of the object should have the longest axis of the object aligned with the x axis, the second longest with the y axis and the shortest with the z axis.

8.2.1 The orientation of pixel detectors

The location and orientation of pixel detectors is particularly important for diffraction experiments.

For convenience specific rules have been defined for the coordinate system of pixel detectors. Instead of defining a local origin plus a default orientation the detectors are defined in absolute terms by a `corner_position` plus a set of `basis_vectors`.

The `corner_position` should contain the x, y and z coordinates of the corner of the first data element which corresponds to the corner of the detector. The `corner_position` must always be defined.

Figure 11 gives an example.




Figure 11: The coordinates of the four corners of the first pixel in a pixel detector.

The `basis_vectors` are a matrix containing a set of 3D vectors from the center of the first element to the center of the second element for each dimension of the detector. The number of rows of the matrix (first dimension) is equal to the number of dimensions of data and the number of columns (second dimension) is equal to 3 (the x, y and z position in 3D space).

The first row then defines the relative distance between the first two pixels in the first dimension (e.g. position of (1,0) - position of (0,0) for a 2D array), and the second row elements do the same along the second dimension (e.g. position of (0,1) - position of (0,0) for a 2D array).

If no `basis_vectors` are specified they are assumed to be:

$$\begin{bmatrix} 0 & -y_pixel_size & 0 \\ -x_pixel_size & 0 & 0 \end{bmatrix}$$

This results in the first dimension parallel to the y axis and the second to the x axis. This convention should be used as much as possible as many programs will assume it when displaying data. If more dimensions are required it's strongly encouraged to make the last two dimension correspond to the y and x axis.




Figure 12: Pixel locations for a pixel detector with no basis_vectors defined or rotation applied and with a corner_position of (pixel width, pixel height, detector distance).

9 Scans

While for simplicity it preferable to keep only one image in each Entry class, often there are experiments where many images come together to form a single dataset. A common case is when one is obtaining multiple images of the same sample at different rotations to performance a tomographic reconstruction. We designate datasets such as these, where one or more variables are changed while images are collected, as *scans*.

In a scan the multiple images are stored in a single 3D or higher dimensional dataset. This has both space and performance advantages, besides making sure the data is kept together. The variables being scanned should always be the first dimensions with the last dimensions reserved for the physical dimensions of the detector. For example in the case of a ptychographic dataset, where the transla-

tion of the sample is scanned, the dimensions of the data would be translation, y dimension of the CCD, x dimension of the CCD.

9.1 Dimension Scales

When using scans it is crucial to determine what variable corresponds to each dimension. To define this you need to use the *axes attribute*. The `axes` attribute should specify a colon separated list of names of the datasets that correspond to each dimension. The named dataset, or a link to it, should exist in the same group as the field with the `axes` attribute. The exception to this rule are the names which match the *implicit axes* of the dataset. For example “y” and “x” which are assumed to correspond to the number of pixels of the detector in each dimension, for many datasets. For details about the *implicit axes* see the [next subsection](#).

For example in the case of tomography the `axes` attribute of the `data` field could have the value `orientation:y:x` and in the same group there should be a link named `orientation` that links to the `orientation` field of the sample being imaged. The “x” and “y” are implicitly assumed to correspond to the detector dimensions.

In another important use case, storing large XFEL datasets, the `axes` attribute of the `data` field could have the value `experiment_identifier:y:x`, and a link named `experiment_identifier` pointing to the `experiment_identifier` field at the [Entry](#) level. The `experiment_identifier` would then be a list of strings, one per shot, with a unique identifier for each of the shots.

9.2 Implicit Axes

Every dataset with multiple dimensions has a set of *implicit axes* which describe what each dimension of the dataset correspond to. These implicit axes are described in the CXI entry reference for the group that contains the dataset. For example the [Detector class](#) contains a `corner_position` dataset with implicit axes = `coordinate`. For datasets which have more dimensions than implicit axes, the extra axes should be specified using the `axes` attribute as described [above](#). For datasets with less dimensions than implicit axes, only the last axes apply.

9.3 Modular Pixel Detectors

Nowadays many pixel array detectors are composed of many small modules tiled together, such as the (in)famous CSPAD at LCLS (see [Fig. 13](#)).

The best way to save such detector images in a CXI file is to populate the “module_identifier” field of the Detector with a string identifier for each of the modules and add an extra dimension to all other fields, using 9 scans, which would identify the module in question.

For example if we wished to save an image from a two piece detector we could set the module_identifier to [“lower-half”, “upper-half”], set the axes attribute of the data field to module_identifier:y:x, and set the axes attribute of the corner_position field to module_identifier:coordinate




Figure 13: The Cornell-SLAC Pixel Area Detector (CSPAD), an example a modular pixel array detector. Note that the four quadrant can move with respect to each other changing the size of the hole in the middle. Thus the detector geometry is variable.

When storing large datasets with many images from modular detectors it’s natural to combine the modular pixel detector feature with an image stack scan to create a 4D array representing the multiple 2D modules that make each image inside of a large stack of images. For an example see Fig. 6.

A CXI entries reference

A.1 Top level (root)

This node represents the top level of the HDF5 file and holds some general information about the file as well as number of entries.

Table 2: CXI top level entries

Member	Type	Quantity
<code>cxi_version</code>	int	version
<code>entry_N</code>	Entry class	
<code>number_of_entries</code>	int	unitless

`cxi_version` - CXI format version times 100. Version 1.00 would be represented by 100 and 1.5 by 150.

`entry_N` - The measurements recoded in this file.

`number_of_entries` - Total number of entries in the file.

A.2 Attenuator

This class describes a beamline attenuator used during data collection.

Table 3: Attenuator class members

Member	Type	Quantity
<code>distance</code>	float	length
<code>thickness</code>	float	length
<code>attenuator_transmission</code>	float	unitless
<code>type</code>	string	text

`distance` - Distance from sample.

`thickness` - Thickness of attenuator along beam direction.

`attenuator_transmission` - The nominal amount of the beam that gets through (transmitted intensity)/(incident intensity).

`type` - Type or composition of attenuator.

A.3 Beam

This class describes the properties of the X-ray beam during data collection.

Table 4: Beam class members

Member	Type	Quantity
incident_energy	float array	variable
incident_wavelength	float array	variable
incident_polarization	float array	variable

[incident_energy](#) - Energy of the incident beam.

[incident_wavelength](#) - Wavelength of the incident beam.

[incident_polarization](#) - Polarization of the incident beam.

A.4 Data

This class is a general placeholder for the most important information in each Entry class. It is mandatory that there is at least one Data class in each Entry class. Most data analysis and plotting programs will primarily focus in this class.

Table 5: Data class members

Member	Type	Quantity
data	float/complex array	variable
errors	float/complex array	variable

[data](#) - Most important data values. Implicit axes = $y : x$.

[errors](#) - Standard deviations of data values. Implicit axes = $y : x$.

A.5 Detector

This class holds information about one of the detectors used during the experiment. Raw data recorded by a detector as well as its position and geometry should be stored in this class.

[basis_vectors](#) - A matrix with the basis vectors of the detector data. For more details see 8.2.1. Implicit axes = `dimension:coordinate`.

[corner_position](#) - The x, y and z coordinates of the corner of the first data element. For more details see 8.2.1. Implicit axes = `coordinate`.

[counts_per_joule](#) - Number of counts recorded per each joule of energy received by the detector. The number of incident photons can then be calculated

Table 6: Detector class members

Member	Type	Quantity
basis_vectors	$N \times 3$ float array	length
corner_position	3 floats	length
counts_per_joule	float	unitless
data	float array	variable
data_dark	float array	variable
data_error	float array	variable
data_sum	float	variable
data_white	float array	variable
description	string	text
distance	float_array	variable
start_time	string	variable
geometry_1	Geometry class	
mask	int array	unitless
module_identifier	string array	text
score	float array	variable
tags	boolean array	unitless
calibration_date	string	date
calibration_method	string	text
x_pixel_size	float	length
y_pixel_size	float	length

by:

$$\text{number of photons} = \frac{\text{data counts}}{\text{source energy} \times \text{counts per joule}}$$

data - Recorded signal values. Implicit axes = $y : x$.

data_dark - Image recorded with the shutter closed, used for background subtractions. Implicit axes = $y : x$.

data_error - The best estimate of the uncertainty in the data value. Where possible, this should be the standard deviation, which has the same units as the data. Implicit axes = $y : x$.

data_sum - Sum of all the elements in the data array. This number is often useful as a cheap measure of data quality.

data_white - Image recorded without the sample, used for background subtractions. Implicit axes = $y : x$.

description - name/manufacturer/model/etc. information.

distance - Closest distance from the detector to the sample for each of the mea-

surements/images.

geometry_1 - Position and orientation of the center of mass of the detector. This should only be specified for non pixel detectors. For pixel detectors use **basis_vectors** and **corner_position**.

mask - Not all the pixels in a detector might have the same value. This 32bit mask makes it possible to distinguish different kinds of pixels.

start_time - Time stamp for each of the measurement/images. The time stamp must follow the ISO8601 format with nanosecond precision.

The following list defines the meaning of each bit when is it set, as well as the names of constants, defined in `cxih.h`, useful for checkings their values:

Table 7: Definition of each bit in the mask

Bit	Meaning	Constant
0x00000001	If set the pixel is invalid	CXI_PIXEL_IS_INVALID
0x00000002	If set the pixel is saturated	CXI_PIXEL_IS_SATURATED
0x00000004	If set the pixel is hot, meaning it typically reads a high value unrelated to the signal.	CXI_PIXEL_IS_HOT
0x00000008	If set the pixel is dead, meaning it typically reads a low value unrelated to the signal.	CXI_PIXEL_IS_DEAD
0x00000010	If set the pixel is under a shadow	CXI_PIXEL_IS_SHADOWED
0x00000080	If set the pixel value should not be trusted, for undefined reasons.	CXI_PIXEL_IS_BAD
0x00000200	If set the pixel does not correspond to a real pixel, but to a gap region in a modular detector.	CXI_PIXEL_IS_MISSING
0x00000400	If set the pixel contains unusually high levels of noise.	CXI_PIXEL_IS_NOISY
0x00001000	If set the pixel signal is above the background	CXI_PIXEL_HAS_SIGNAL

All other bits have no standard meaning and can be used for any purpose the user sees fit. More bits will be defined as the format evolves so users are encouraged to use the high bits to avoid collisions. A pixel with no bits sets correspond to an ideal pixel. Implicit axes = $y : x$.

module_identifier - A string identifier for each of the modules that constitute the pixel detector. See also [Modular Pixel Detector](#).

score - The score, or an array of scores for each of the measurements/images. The fastest dimension should match the number of measurements and

the slowest the number of scores used. Dimensions with only one element can be collapsed (e.g. if only one score type is used the array can be one dimensional). Should have a `headings` attribute containing an array with the names of all the scores. `score` is useful, for example, to give a hit score to each image in a stack.

tags - The tags which apply to each of the measurements/images. The fastest dimension should match the number of measurements and the slowest the number of `tags` used. Dimensions with only one element can be collapsed (e.g. if there is only one image with multiple tags, the `tags` array can be one dimensional). Each entry should be either a 1, signifying that the measurement is tagged with that particular tag, or 0 otherwise. As HDF5 currently does not have native boolean support the dataset should be saved as an 8-bit integer, e.g. `H5T_NATIVE_B8`. It should have a `headings` attribute containing an array with the names of all the `tags`. Tags are useful for organizing and filtering large datasets.

x_pixel_size - Width of each detector pixel.

y_pixel_size - Height of each detector pixel.

A.6 Entry

Base CXI class which holds all other classes.

Table 8: Entry class members

Member	Type	Quantity
<code>data_N</code>	Data class	
<code>end_time</code>	string	text
<code>experiment_identifier</code>	string	text
<code>experiment_description</code>	string	text
<code>image_N</code>	Image class	
<code>instrument_N</code>	Instrument class	
<code>program_name</code>	string	text
<code>sample_N</code>	Sample class	
<code>start_time</code>	string	text
<code>title</code>	string	text

data_N - Main data collected.

end_time - Ending time of measurement.

experiment_identifier - Unique identifier for the experiment, defined by the facility, possibly linked to the proposals, and each image in the dataset.

`experiment_description` - Description of the experiment.

`image_N` - Processed images.

`instrument_N` - Instrument used.

`program_name` - Name of program used to generate this file.

`sample_N` - Sample used.

`start_time` - Starting time of measurement.

`title` - Extended title for entry

A.7 Geometry

This class holds the general position and orientation of a component.

Table 9: Geometry class members

Member	Type	Quantity
<code>orientation</code>	6 floats	unitless
<code>translation</code>	3 floats	length

`orientation` - Dot products between the local and the global unit vectors. Implicit axes = `dot_product`.

`translation` - The x, y and z components of the translation of the origin of the object relative to the origin of the global coordinate system (the place where the X-ray beam meets the sample). Implicit axes = `coordinate`.

Only one orientation and one translation is permitted in each geometry class.

The position of the origin of the object should be explicitly defined for each object. If it is not defined it should be assumed to be the center of the object.

The orientation information is stored as direction cosines. The direction cosines will be between the local coordinate directions and the global coordinate directions. The unit vectors in both the local and global coordinates are right-handed and orthonormal.

Calling the local unit vectors (x', y', z') and the reference unit vectors (x, y, z) the six numbers will be $[x' \cdot x, x' \cdot y, x' \cdot z, y' \cdot x, y' \cdot y, y' \cdot z]$ where " \cdot " is the scalar dot product (cosine of the angle between the unit vectors).

Notice that this corresponds to the first two rows of the rotation matrix that transforms from the global orientation to the local orientation, namely, the matrix that when applied to x' produces the vector $(1, 0, 0)$ and analogously for y'

and z' . It also corresponds to the first two *columns* of the matrix that when applied to each of the global unit vector (x, y, z) produces respectively (x', y', z') . The third row can be recovered by using the fact that the basis vectors are orthonormal.

A.8 Image

This class should be used to store processed image data. It describes what analysis has been done, as well as holding important information for further image processing. It should not be used for raw data storage, which should be stored in the Detector class.

Table 10: Image class members

Member	Type	Quantity
title	string	text
data	float/complex array	variable
data_error	float/complex array	variable
data_space	string	text
data_type	string	text
detector_N	Detector class	
dimensionality	int	unitless
image_center	3 floats	pixels
image_size	3 floats	(inverse) length
is_fft_shifted	int	unitless
mask	int array	unitless
process_N	Process class	
reciprocal_coordinates	3 floats array	inverse length
source_N	Source class	

title - A title for the image.

data - The value of the image at each pixel. Implicit axes = $y : x$.

data_error - The best estimate of the uncertainty in the data value. Where possible, this should be the standard deviation, which has the same units as the data. Implicit axes = $y : x$.

data_space - Specifies if the image lives in real or diffraction (Fourier) space. Only has two valid values: “real” and “diffraction”.

data_type - Defines what the data represents. The following values are allowed: “intensity”, “electron density”, “amplitude”, “unphased amplitude”, “autocorrelation”. “amplitude” implies a phased dataset, while “unphased amplitude” corresponds to the square root of the intensity.

`detector_N` - Link to the detectors used to obtain this image.

`dimensionality` - Number of dimensions of the image. Restricted to 1 2 or 3.

`image_center` - The location of the zero frequency component on a diffraction image in fractional pixels (see Fig. 14 for the coordinate system convention).




Figure 14: Pixel coordinate system for a 2x2 pixels image.

Implicit axes = `coordinate`.

`image_size` - The width, height and depth of the image. For real space images this corresponds to the length and for reciprocal space ones to the inverse length of the sides of the image. Implicit axes = `dimension`.

`is_fft_shifted` - If set to 1 the image is assumed to have the quadrants shifted (see Fig. 15).




Figure 15: FFT shifting a 2 dimensional image.

`mask` - 32-bit unsigned integer mask specifying the properties of each pixel.

The following list defines the meaning of each bit when is it set, as well as the names of constants, defined in `cxih`, useful for checkings their values:

All other bits have no standard meaning and can be used for any purpose

Table 11: Definition of each bit in the mask

Bit	Meaning	Constant
0x00000001	If set the pixel is invalid	CXI_PIXEL_IS_INVALID
0x00000002	If set the pixel is saturated	CXI_PIXEL_IS_SATURATED
0x00000004	If set the pixel is hot	CXI_PIXEL_IS_HOT
0x00000008	If set the pixel is dead	CXI_PIXEL_IS_DEAD
0x00000010	If set the pixel is under a shadow	CXI_PIXEL_IS_SHADOWED
0x00001000	If set the pixel signal is above the background	CXI_PIXEL_HAS_SIGNAL
0x00010000	If set the pixel is inside of the reconstruction support	CXI_PIXEL_INSIDE_SUPPORT

the user sees fit. More bits will be defined as the format evolves so users are encouraged to use the high bits to avoid collisions.

Implicit axes = $y : x$

process_N - Processes used to obtain this image. They should be listed in chronological order with the first processed used named process₁, the second process₂ and so on.

reciprocal_coordinates - The diffraction (Fourier) space coordinates of the center of each pixel. Note that the dimension corresponding to the 3 different components should go before the image dimensions. So for example for an image of size [10,20,5] te reciprocal_coordinates will have size [3,10,20,5].

Implicit axes = coordinate : $y : x$

source_N - Link to the source used to obtain this image.

A.9 Instrument

Template of instrument descriptions comprising various beamline components. Each component will also be a class defined by its distance from the sample. Negative distances represent beamline components that are before the sample while positive distances represent components that are after the sample. This device allows the unique identification of beamline components in a way that is valid for both reactor and pulsed instrumentation.

Each Instrument instance corresponds to one beamline.

name - Name of the instrument.

attenuator_N - The attenuators that are part of the instrument.

detector_N - The detectors that compose the instrument.

Table 12: Instrument class members

Member	Type	Quantity
name	string	text
attenuator_ <i>N</i>	Attenuator class	
detector_ <i>N</i>	Detector class	
monochromator_ <i>N</i>	Monochromator class	
source_ <i>N</i>	Source class	

[source_*N*](#) - The source used by the instrument.

A.10 Monochromator

Define a monochromator used in the instrument.

Table 13: Monochromator class members

Member	Type	Quantity
energy	float	energy
energy_error	float	energy

[energy](#) - Peak of the spectrum that the monochromator selects.

[energy_error](#) - Standard deviation of the spectrum that the monochromator selects.

A.11 Log

This class can be used to store time-stamped information recorded during a run. For example, it can be used to store readouts from environmental sensors (temperature, pressure) or event codes from a triggering system.

Table 14: Log class members

Member	Type	Quantity
time	char array	variable
value	char/float/int array	variable
raw_value	char/float/int array	variable
description	string	text

[time](#) - Time stamp of the logged entry. Must follow the ISO 8601 format with nanosecond precision.

value - Logged entry.

raw_value - Raw sensor readout used to compute the logged entry. For example thermocouple voltage of a temperature sensor.

description - Description of the logged value.

A.12 Note

This class can be used to store additional information in a CXI file e.g. additional text logs, configuration files, pictures, movies, audio.

Table 15: Note class members

Member	Type	Quantity
author	string	text
data	char array	binary/text
date	string	text
description	string	text
file_name	string	text
type	string	text

author - Author or creator of note.

data - Binary/text note data.

date - Date note created/added.

description - Title of an image or other details of the note

file_name - Name of original file name if note was read from an external source.

type - Mime content type of note data field e.g. image/jpeg, text/plain, text/html.

A.13 Process

Document an event of data processing, reconstruction, or analysis.

command - Command line used to run the program.

comments - Comments related to how the data was processed.

date - Date and time of processing in ISO 8601 format.

note_N - Notes providing extra information like configuration files used, other inputs required or any other important information.

Table 16: Process class members

Member	Type	Quantity
command	string	text
comments	string	text
date	string	text
note_ <i>N</i>	Note class	
program	string	text
version	string	text

program - Name of the program used.

version - Version of the program used.

A.14 Result

This class should be used to store the results of an analysis procedure which produces non-image data (image data should be stored in the more specific [Image class](#)). It describes what analysis has been done, as well as holding important information for further processing. It should not be used for raw data storage, which should be stored in the Detector class.

Table 17: Result class members

Member	Type	Quantity
data	float/complex array	variable
data_error	float/complex array	variable
description	string	text
detector_ <i>N</i>	Detector class	
experiment_identifier	string	text
nPeaks	int array	unitless
peakMaximumValue	float array	variable
peakNPixels	int array	unitless
peakSNR	float array	variable
peakTotalIntensity	float array	variable
peakXPosAssembled	float array	variable
peakYPosAssembled	float array	variable
peakXPosRaw	float array	variable
peakYPosRaw	float array	variable
data	float/complex array	variable
process_ <i>N</i>	Process class	
source_ <i>N</i>	Source class	

data - The value of the image at each pixel. Implicit axes = $y : x$.

data_error - The best estimate of the uncertainty in the data value. Where possible, this should be the standard deviation, which has the same units as the data. Implicit axes = $y : x$.

description - Short description of the analysis performed.

detector_N - Link to the detectors used to obtain the input data for this analysis.

experiment_identifier - Unique identifier for the experiment, defined by the facility, possibly linked to the proposals, and each image in the dataset. Usually a symlink to the one in the respective [Entry class](#)

nPeaks - Number of peaks in found in each image. 1D array with a length equal to the number of images. Implicit axes = `experiment_identifier`.

peakMaximumValue - The maximum value within the peak. 2D array with the first (slow varying) dimension matching the number of images and the second (fast varying) at least as big as the number of peaks in the image. Any extra values should be ignored. Implicit axes = `experiment_identifier:nPeaks`.

peakNPixels - The number of pixels within the peak. 2D array with the first (slow varying) dimension matching the number of images and the second (fast varying) at least as big as the number of peaks in the image. Any extra values should be ignored. Implicit axes = `experiment_identifier:nPeaks`.

peakSNR - The signal to noise ratio in the peak. 2D array with the first (slow varying) dimension matching the number of images and the second (fast varying) at least as big as the number of peaks in the image. Any extra values should be ignored. Implicit axes = `experiment_identifier:nPeaks`.

peakTotalIntensity - The total intensity in the peak. 2D array with the first (slow varying) dimension matching the number of images and the second (fast varying) at least as big as the number of peaks in the image. Any extra values should be ignored. Implicit axes = `experiment_identifier:nPeaks`.

peakXPosAssembled - The x position of the peak in the assembled image. 2D array with the first (slow varying) dimension matching the number of images and the second (fast varying) at least as big as the number of peaks in the image. Any extra values should be ignored. Implicit axes = `experiment_identifier:nPeaks`.

peakYPosAssembled - The y position of the peak in the assembled image. 2D array with the first (slow varying) dimension matching the number of images and the second (fast varying) at least as big as the number of peaks in the image. Any

extra values should be ignored. Implicit axes = `experiment_identifier:nPeaks`.

`peakXPosRaw` - The x position of the peak in the raw image. 2D array with the first (slow varying) dimension matching the number of images and the second (fast varying) at least as big as the number of peaks in the image. Any extra values should be ignored. Implicit axes = `experiment_identifier:nPeaks`.

`peakYPosRaw` - The y position of the peak in the raw image. 2D array with the first (slow varying) dimension matching the number of images and the second (fast varying) at least as big as the number of peaks in the image. Any extra values should be ignored. Implicit axes = `experiment_identifier:nPeaks`.

`process_N` - Processes used to obtain this result. They should be listed in chronological order with the first processed used named `process_1`, the second `process_2` and so on.

`source_N` - Link to the source used to obtain the input data for this analysis.

A.15 Sample

This class holds basic information about the kind of sample used, its geometry and properties.

Table 18: Sample class members

Member	Type	Quantity
<code>concentration</code>	float	mass/volume
<code>description</code>	string	text
<code>geometry_1</code>	Geometry class	
<code>mass</code>	float	mass
<code>name</code>	string	text
<code>temperature</code>	float	temperature
<code>unit_cell</code>	6 floats	length/angle
<code>unit_cell_group</code>	string	text
<code>thickness</code>	float	length
<code>unit_cell_volume</code>	float	volume

`concentration` - Concentration of the sample.

`description` - Description of the sample.

`geometry_1` - Position and orientation of the center of mass of the sample.

`mass` - Mass of sample.

`name` - Descriptive name of sample.

`temperature` - Sample temperature.

`unit_cell` - Unit cell parameters (a,b,c α, β, γ). Implicit axes = `unit_cell`.

`unit_cell_group` - Crystallographic space group of the crystal in PDB format.

`thickness` - Sample thickness.

`unit_cell_volume` - Volume of the unit cell.

A.16 Source

Class describing the light source being used.

Table 19: Source class members

Member	Type	Quantity
<code>energy</code>	float	energy
<code>name</code>	string	text
<code>pulse_energy</code>	float	energy
<code>pulse_width</code>	float	time

`energy` - Energy of each photon.

`name` - The name of the source, for example ALS.

`pulse_energy` - Sum of the energy of all the photons in the pulse.

`pulse_width` - Duration of the pulse.

B Changes from previous versions

B.1 Version 1.5

- Updated the mask bits (see [Image class](#) and [Detector class](#)) to match common usage (in cheetah and owl).
- Added peak lists to the [Result class](#), for marking the location of reflections.

B.2 Version 1.4

- Introduced the concept of [implicit axes](#).
- Added support for [modular pixel detectors](#).
- Added the `module_identifier` field to the [Detector class](#) to support modular pixel detectors.
- Added a [Result class](#) to store non-image analysis results.
- Added more examples, showing the use of scans and modular detectors.
- Moved the examples in the appendix to the main text.
- Added the `tags` and `score` fields [Detector class](#) to aid with filtering and sorting large dataset.

B.3 Version 1.3

- Introduced [Scans](#).
- Added the `axes` attribute to describe dataset dimensions during scans.
- Added the fields `data_dark` and `data_white` to the [Detector class](#).
- Made the type of all the fields that start with `data` in the [Detector class](#) float instead of float/complex.
- Merged the Translation and Orientation classes in the [Geometry class](#).
- Correct formula for number of photons in [Detector class](#).
- Updated the examples to correspond to the changes, in particular the ptychography example now makes use if scans.
- Added a `name` field to the [Source class](#).
- Clarify the description of [Geometry orientation](#).

B.4 Version 1.2.1

- Added a [Monochromator](#) class.

B.5 Version 1.2

- A new entry has been added to the mask of the [Image](#) class to represent the support used during the reconstruction of an image.

B.6 Version 1.1

- Add `pulse_energy` to the [Source](#) class.
- Add `data_sum` and change `counts_per_eV` to `counts_per_joule` in the [Detector](#) class.
- Added a section (this section) listing all the changes from the previous versions.

B.7 Version 1.0

- All quantities must now be in SI units, without exceptions.
- The [geometry of CCD detectors](#) is now represented by a series of basis vectors representing the sides of the detector, as well as the absolute corner position, instead of rotations and translations.

C Diagrams color code




Figure 16: Explanation of the color code used in the diagrams

D Code examples

All the code examples as well as the resulting CXI files are available from cxidb.org.

D.1 Creating a minimal CXI file

Listing 2: Creating a minimal CXI file

```
#!/usr/bin/env python

import h5py
import numpy as np
import math

fileName = "minimal.cxi"
# open the HDF5 CXI file for writing
f = h5py.File(fileName, "w")

# create data
x = np.arange(-5, 5, 0.1)
y = np.arange(0, 5, 0.1)
xx, yy = np.meshgrid(x, y)
sinc = np.sin(xx**2+yy**2)/(xx**2+yy**2)

# populate the file with the classes tree
entry_1 = f.create_group("entry_1")
data_1 = entry_1.create_group("data_1")
# write the data
data = data_1.create_dataset("data", data=sinc)

f.close()
```

The resulting file should be equivalent to the one in Fig. 1.

D.2 Creating a typical raw CXI file

Listing 3: Creating a typical raw CXI file

```
#!/usr/bin/env python

import h5py
import numpy as np
import math

fileName = "typical_raw.cxi"
# open the HDF5 CXI file for writing
f = h5py.File(fileName, "w")
f.create_dataset("cxi_version", data=120)

# create data 1
x = np.arange(-5, 5, 0.1)
y = np.arange(-5, 5, 0.2)
xx, yy = np.meshgrid(x, y)
sinc1 = np.sin(xx**2+yy**2)/(xx**2+yy**2)

# create data 2
x = np.arange(-1, 1, 0.02)
y = np.arange(-1, 1, 0.04)
xx, yy = np.meshgrid(x, y)
sinc2 = np.sin(xx**2+yy**2)/(xx**2+yy**2)

# populate the file with the classes tree
entry_1 = f.create_group("entry_1")
entry_1.create_dataset("experimental_identifier", data=
    "LCLS_2009_Dec11_170451_21963")
entry_1.create_dataset("start_time", data=
    "2009-12-11T17:04:51-0800")
sample_1 = entry_1.create_group("sample_1")
sample_1.create_dataset("name", data="Mimivirus")
instrument_1 = entry_1.create_group("instrument_1")
instrument_1.create_dataset("name", data="AMO")
source_1 = instrument_1.create_group("source_1")
source_1.create_dataset("energy",
    data=2.8893e-16) # in J
source_1.create_dataset("pulse_width",
    data=70e-15) # in s

detector_1 = instrument_1.create_group("detector_1")
detector_1.create_dataset("distance",
    data=0.15) # in meters
detector_1.create_dataset("data", data=sinc1)
```

```
detector_2 = instrument_1.create_group("detector_2")
detector_2.create_dataset("distance",
                          data=0.65) # in meters
detector_2.create_dataset("data", data=sinc2)

data_1 = entry_1.create_group("data_1")
data_1["data"] = h5py.SoftLink('/entry_1/instrument_1/
                               detector_1/data')

data_2 = entry_1.create_group("data_2")
data_2["data"] = h5py.SoftLink('/entry_1/instrument_1/
                               detector_2/data')

f.close()
```

The resulting file should be equivalent to the one in Fig. 2.

D.3 NeXus compatible version of the typical raw CXI file

Listing 4: NeXus compatible version of the typical raw CXI file

```
#!/usr/bin/env python

import h5py
import numpy as np
import math

fileName = "nexus.cxi"
# open the HDF5 CXI file for writing
f = h5py.File(fileName, "w")

f.attrs['NeXus_version'] = "4.3.0"
f.attrs['creator'] = "CXI2NeXus"

f.create_dataset("cxi_version", data=120)

# create data 1
x = np.arange(-5, 5, 0.1)
y = np.arange(0, 5, 0.1)
xx, yy = np.meshgrid(x, y)
sinc1 = np.sin(xx**2+yy**2)/(xx**2+yy**2)

# create data 2
x = np.arange(-1, 1, 0.02)
y = np.arange(0, 1, 0.02)
xx, yy = np.meshgrid(x, y)
sinc2 = np.sin(xx**2+yy**2)/(xx**2+yy**2)

# populate the file with the classes tree
entry_1 = f.create_group("entry_1")
entry_1.create_dataset("experimental_identifier", data=
    "LCLS_2009_Dec11_170451_21963")
entry_1.create_dataset("start_time", data=
    "2009-12-11T17:04:51-0800")
entry_1.attrs['NX_class'] = "NXentry"

sample_1 = entry_1.create_group("sample_1")
sample_1.create_dataset("name", data="Mimivirus")
sample_1.attrs['NX_class'] = "NXsample"

instrument_1 = entry_1.create_group("instrument_1")
instrument_1.create_dataset("name", data="AMO")
instrument_1.attrs['NX_class'] = "NXinstrument"

source_1 = instrument_1.create_group("source_1")
energy = source_1.create_dataset("energy", data=2.8893e-16) # in
```

```
J
energy.attrs['units'] = "J"
pulse_width = source_1.create_dataset("pulse_width",
                                     data=70e-15) # in s

pulse_width.attrs['units'] = "s"
source_1.attrs['NX_class'] = "NXsource"

detector_1 = instrument_1.create_group("detector_1")
distance = detector_1.create_dataset("distance",
                                    data=0.15) # in meters

distance.attrs['units'] = "m"
data = detector_1.create_dataset("data", data=sinc1)
data.attrs['signal'] = 1
data.attrs['units'] = "counts"
detector_1.attrs['NX_class'] = "NXdetector"

detector_2 = instrument_1.create_group("detector_2")
distance = detector_2.create_dataset("distance", data=0.65) # in
    meters
distance.attrs['units'] = "m"
data = detector_2.create_dataset("data", data=sinc2)
data.attrs['units'] = "counts"
detector_2.attrs['NX_class'] = "NXdetector"

data_1 = entry_1.create_group("data_1")
data_1["data"] = h5py.SoftLink('/entry_1/instrument_1/
    detector_1/data')
data_1.attrs['NX_class'] = "NXdata"

data_2 = entry_1.create_group("data_2")
data_2["data"] = h5py.SoftLink('/entry_1/instrument_1/
    detector_2/data')
data_2.attrs['NX_class'] = "NXdata"

f.close()
```

The resulting file should be equivalent to the one in Fig. 9.

D.4 3D complex valued Image file

Listing 5: 3D complex valued Image file

```
#!/usr/bin/env python

import h5py
import numpy as np
import math

fileName = "phased_3d.cxi"
# open the HDF5 CXI file for writing
f = h5py.File(fileName, "w")
f.create_dataset("cxi_version", data=120)

# create data
zz, yy, xx = np.mgrid[-5:5:8j, -5:5:12j, -5:5:16j]
sinc = np.sin(xx**2+yy**2+zz**2)/(xx**2+yy**2+zz**2) + \
        1j*np.cos(xx**2+yy**2+zz**2)/(xx**2+yy**2+zz**2)

# populate the file with the classes tree
entry_1 = f.create_group("entry_1")
sample_1 = entry_1.create_group("sample_1")
sample_1.create_dataset("name", data="Mimivirus")
image_1 = entry_1.create_group("image_1")
image_1.create_dataset("data", data=sinc)
image_1.create_dataset("data_type", data="electron density")
image_1.create_dataset("data_space", data="real")
image_1.create_dataset("image_size",
                      data=[1.65e-6, 1.65e-6, 1.65e-6])
source_1 = image_1.create_group("source_1")
source_1.create_dataset("energy", data=2.8893e-16) # in J
detector_1 = image_1.create_group("detector_1")
detector_1.create_dataset("distance",
                        data=0.15) # in meters
detector_1.create_dataset("x_pixel_size",
                        data=15e-6) # in meters
detector_1.create_dataset("y_pixel_size",
                        data=15e-6) # in meters
process_1 = image_1.create_group("process_1")
process_1.create_dataset("command", data="find_center -i"
                        " mimi_raw.cxi -o mimi_center.cxi")
process_2 = image_1.create_group("process_2")
process_2.create_dataset("command", data="phase_image"
                        " phasing.txt -i mimi_center.cxi"
                        " -o mimi_phased.cxi")
note_1 = process_2.create_group("note_1")
note_1.create_dataset("file_name", data="phasing.txt")
```

```
note_1.create_dataset("description", data="configuration text"
                      " file used for phasing")
note_1.create_dataset("data", data='algorithm = "HIO"')

data_1 = entry_1.create_group("data_1")
data_1["data"] = h5py.SoftLink('/entry_1/image_1/data')

f.close()
```

The resulting file should be equivalent to the one in Fig. 8.

D.5 Spartan viewer for CXI files

Listing 6: Spartan viewer for CXI files

```
#!/usr/bin/python
import h5py
import sys
from PyQt4 import QtGui, QtCore, Qt
import matplotlib.pyplot as plt
from operator import mul
import numpy
import signal
import scipy.interpolate
import scipy.ndimage

signal.signal(signal.SIGINT, signal.SIG_DFL)

def onclick(event):
    z = plt.gca().get_images()[0].get_array()
    x = numpy.arange(0, z.shape[0], 1)
    y = numpy.arange(0, z.shape[1], 1)
    v = scipy.ndimage.map_coordinates(z, [[event.ydata], [event
        .xdata]], order=1)
    print 'xdata=%f, ydata=%f value=%e'%(
        event.xdata, event.ydata, v[0])

class Viewer(QtGui.QMainWindow):
    def __init__(self):
        QtGui.QMainWindow.__init__(self)
        self.tree = QtGui.QTreeWidget(self)
        self.setCentralWidget(self.tree)
        self.buildTree()
        self.tree.itemClicked.connect(self.handleClick)
    def handleClick(self, item, column):
        if(item.text(column) == "Click to display"):
            data = self.datasets[str(item.text(2))]
            fig = plt.figure()
            ax = fig.add_axes([0, 0, 1, 1])
            if(numpy.iscomplexobj(data)):
                data = numpy.abs(data)
            if(len(data.shape) == 1):
                plt.plot(data)
            else:
                ax.imshow(data)
            cid = fig.canvas.mpl_connect('
                button_press_event', onclick)
    def buildTree(self):
        self.datasets = {}
        self.tree.setColumnCount(2)
```

```
self.f = h5py.File(sys.argv[1], "r")
item = QtGui.QTreeWidgetItem(QtCore.QStringList("/"))
self.tree.addTopLevelItem(item)
self.buildBranch(self.f, item)
def buildBranch(self, group, item):
    for g in group.keys():
        lst = QtCore.QStringList(g)
        if(isinstance(group[g], h5py.Group)):
            child = QtGui.QTreeWidgetItem(lst)
            self.buildBranch(group[g], child)
            item.addChild(child)
        else:
            if(not group[g].shape or reduce(mul, group[g]
                ].shape) < 10):
                lst.append(str(group[g][()]))
            else:
                lst.append("Click to display")
                lst.append(group[g].name)
                self.datasets[group[g].name] = group[g]
            item.addChild(QtGui.QTreeWidgetItem(lst))

app = QtGui.QApplication(sys.argv)
aw = Viewer()
aw.show()
sys.exit(app.exec_())
```

This simple program provides a way to view CXI files. Does not support 3D or complex valued data. To use it pass the name of the file you which to view as a command line argument.