





## How to take a whooping cough oral fluid swab







Write the date of sample on the label (G) and on the request form (E). Stick the label to the clear tube (B)







Place clear tube (B) into self adhesive clear bag (C). Push air out of bag. Remove red strip from top fold down corner and stick down. Now wash your hands.



Place form (E) with clear bag containing the sample (C) into cardboard box (D)



Place cardboard box (D) into pre-paid grey plastic envelope (F) and seal. Post pre-paid grey plastic envelope (F) in your local Royal Mail post box